Gebrauchsanweisung | *User instruction*
actric MUSICearmold

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EN actric MUSICearmold: The customised hearing protection is reusable. It is worn in the ear canal and protects against dangerous or harmful noise when used correctly, i.e. when correctly fitted in the ear. Any other use is not intended and therefore not permitted.

The correct use: Hold the edge or handle with your thumb and forefinger and push the hearing protection into the ear canal with light pressure until it fits comfortably and completely seals the ear canal.

IMPORTANT: The earmold marked 'R' fits in the right ear, the one marked 'L' in the left ear.

Notes on use: The hearing protection only provides its full protective effect if it fits correctly in the ear and is worn continuously in noisy situations. Ensure that the hearing protection is used correctly in accordance with the ma-

nufacturer's instructions and checked regularly (at least every 6 months) for damage and correct fit and function in the ear. Damage to the product or changes to the ear canal can have a negative impact on the protective function. In this case, the earmold should be replaced. If the situation is unclear or if you suspect poor performance, the hearing protection must be checked by qualified personnel. If these recommendations are not followed, the protective effect may be seriously impaired. The product is not subject to a limited shelf life. The lifetime depends on the use (wearing frequency, wearing time) and care of the product and can be up to 5 years.

WARNING: These hearing protection earmolds should not be worn with a cord if there is a risk of the connecting cord becoming caught or torn during wearing. We then recommend removing the cord. Do not remove the hearing protectors from the ear by pulling on the cord, but by gripping the earmold. Sudden or rapid removal of the hearing protectors from the ear canal can damage the eardrum.

Storage and care: Store your hearing protection in the supplied case at a maximum temperature of +55°C. Care for the hearing protection after each use by removing any ear wax residue with the wax stick provided and cleaning the earmolds with water, some mild soap or a damp cloth.

IMPORTANT: The acoustic filter must not come into contact with moisture. It should be removed from the earmold before cleaning with water. Make sure that the hearing protection is always worn with the attenuating filter in the ear!

This product can be negatively affected by certain chemical substances. For further information, please contact the manufacturer.

EU Type Examination: The customised hearing protection actric MUSICearmold has been tested in accordance with the EN 352-2:2020 standard and is certified and controlled in accordance with the PPE Regulation (EU)2016/425 by PZT GmbH, Bismarckstraße 264 B, 26389 Wilhelmshaven, Germany. You can find the declaration of conformity for download at: link.bachmaier.de/ke-actric

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actric GmbH
bachmaier GmbH